

# Mobile Health Revolution: The MobiCare Approach to Patient Support



كلية الهندسة بنها

FACULTY OF ENGINEERING- BENHA



Supervisor:

Dr. Ayman Soliman

كلية الهندسة بنها  
FACULTY OF ENGINEERING- BENHA



# Presentation Content




## Introduction

- Introduction Our Team
- introduction to the project
- Why this idea was chosen (the problem it solves)



## Problem Statement

- What challenges do people with disabilities face with traditional wheelchairs?
- The importance of a smart, remotely controlled and health-monitored wheelchair



## Objectives

- Facilitate mobility and control
- Monitor health conditions




## System Overview

- How the System Works
- Devices used
- Block Diagram



# Presentation Content



## Embedded Section

- Embedded Components
- Problems and Solutions



## Embedded & AI Integration


- Integration between Embedded Systems and Artificial Intelligence

5

6

7

8



## AI Section

- Voice (How to build the model)
- Testing for voice model
- Face direction Tracking
- EYE direction Tracking



## Mobile App Section

- UI
- Technology (Flutter)
- Database (Firebase)



# Presentation Content



## Mobile & Embedded Link

- Integration between Mobile Applications and Embedded Systems



## Conclusion

- Summary of the idea
- Key outcomes and achievements
- Is the project scalable or ready for development?

9

10

11




## Testing & Results

- Photos and videos of the testing phase
- Data samples and system performance evaluation




# Introduction to the Project

In our project, we sought to solve a problem that affects many people: how to help patients and people with special needs move easily without relying on anyone else. We designed and implemented a smart wheelchair, called MobiCare, that uses modern technologies such as artificial intelligence and mobile applications to control the wheelchair's movement in various ways, such as voice, eye tracking, or even a mobile app. The project's primary goal is to improve patients' quality of life and facilitate their movement in a safe and simple manner.



# Why This Idea Was Chosen (Problem It Solves)

Many individuals with physical disabilities face daily challenges in using traditional wheelchairs, especially those who cannot operate them manually. Most available solutions are either too expensive, lack intelligent control, or do not support health monitoring features. This project aims to address these issues by creating a cost-effective, intelligent, and accessible system that improves mobility, enables remote control, and allows caregivers or doctors to monitor the user's health in real-time.






# What Challenges Do People with Disabilities Face with Traditional Wheelchairs?

Lack of physical ability to control the wheelchair manually

No remote-control options

No health monitoring features

There is no way to track the user's vital signs, making it difficult to detect health issues early.



**Lack of safety and independence**

People using traditional wheelchairs often feel dependent on others to help them move, which takes away their sense of independence and privacy in daily life.

**High cost of existing smart solutions**

Advanced smart wheelchairs on the market are often expensive and not accessible for everyone.



# The Importance of a Smart Wheelchair

## Enables remote control through mobile application


Users or caregivers can move the wheelchair without physical effort, enhancing convenience and safety.

## Integrates health monitoring sensors

Vital signs like heart rate or temperature can be tracked in real-time, allowing quick response from caregivers or doctors.

## Improves user independence

Smart features reduce reliance on others, helping users regain a sense of control over their mobility.







# Objectives Of The Project



## Facilitate Mobility and Control

- Remote control via mobile application.
- Support for users with limited motor abilities.
- More flexibility and freedom.





## Monitor Health Conditions

- Real-time health monitoring using sensors.
- Automatic data collection and storage.
- Improved safety and medical response.




# How the System Works ?





**1) Mode Selection:**  
The user selects the operating mode using push buttons connected to the STM32 microcontroller.




## 2) Manual Mode (Joystick):

- When the joystick mode is selected, the STM32 reads the joystick analog values (X, Y) through the ADC.
- Then, it determines the direction based on the joystick position
- After that, The STM32 receives the movement commands and controls the DC motors through the BTS7960 motor driver to move the wheelchair accordingly.

### 3) AI System (Raspberry Pi - Voice/Eye Tracking):


In AI modes, the Raspberry Pi processes voice or eye-tracking data and sends movement commands to the STM32 via UART.



#### 4) Mobile Application (Mobile Control Mode):

In mobile mode, the smartphone app communicates with the ESP8266 over Wi-Fi.

The app sends movement commands to the ESP8266, which forwards them to the STM32 via UART.



## 5) Health Monitoring (ESP8266 + Firebase):

Health sensors (such as heartbeat and temperature) are connected to the ESP8266.



# Devices used

STM32 MCU



DC MOTOR



Joystick Module



Raspberry Pi 4



Raspberry CAM



ESP WIFI Module



BTS7960  
MOTOR Driver



Temp MIx90614




Heartbeat  
Max30102



Microphone




# Block Diagram





# Block Diagram




# STM32 Microcontroller

## Overview:

STM32 is a fast and efficient microcontroller, widely used in embedded systems.

## In Our Project:



- It works as the main brain of the system.
- Reads input signals from the push buttons to detect the user's chosen mode and movement direction.
- Sends control signals to the BTS7960 motor driver, managing the wheelchair's movement in real-time.



# BTS7960 Motor Driver

## Why We Used It?

- Used to control the direction of the DC motors in the wheelchair.
- Contains an H-Bridge circuit that works like a set of electronic switches (transistors) to reverse the motor's direction.
- Can handle high current, making it suitable for powerful wheelchair motors 43A.





# JOYSTICK (Manual Mode)

## How We Read Direction ?

- The joystick outputs voltage values depending on its movement:
- X and Y values range from ~0 to 4095 (if STM32 ADC is 12-bit).
- THRESHOLD\_LOW = 1000
- THRESHOLD\_HIGH = 3000

Direction	ADC Condition
LEFT	X < THRESHOLD_LOW
RIGHT	X > THRESHOLD_HIGH
UP	Y < THRESHOLD_LOW
DOWN	Y > THRESHOLD_HIGH
CENTER	Values in between (no movement)





# Problems and Solutions

# Problems and Solutions




ATmega32A had unstable ADC & limited UART.



STM32F103 has 3 UARTs & better ADC handling (12 bit vs 10 bit).



STM32 runs faster (72MHz vs 16MHz).




Easier to add new features (e.g., sensors, communication) cause of its has more I/O pins.



Cheaper & more available in the market.

# Block Diagram



# Block Diagram





UI

UI

Code

UI Code

Database





UI

# Why Is UI Important?

Planning Ahead

Clear Design Vision



**UI**

# Why Is UI Important?

**Faster Build Time**

**Simplified Development**

# UI

# Why Is UI Important?




## Understanding User Needs

# Overview

UI

# Light Mode Design Screens

# Light Mode Design Screens



## Sign in Screen

### Modern & Accessible Design:

Clean, user-friendly layout with soft colors (white, teal, blue) reflecting trust and innovation.

### Structured Layout:

Central login card with rounded corners and soft shadow. Background features dynamic oval shapes for visual balance.

### Clarity & Simplicity:

- Logo and app name at top for branding.
- Input fields with icons and placeholders.
- Highlighted “Sign in” button in gradient green.
- Subtle links: “Forgot password?”, “Remember me”.


### Visual Hierarchy:

- Bold typography for headings; light font for guidance. Clear button contrast and logical arrangement.

### Accessibility Focus:

- Large, easy-to-tap fields. Icons aid recognition. Strong text/background contrast for readability.


# Light Mode Design Screens



This screen shares many of the same design principles as the Sign-In Screen, such as clarity, consistency, and accessibility. However, it is uniquely tailored to support new user registration, offering additional input fields and features specific to account creation.

# Light

# Mode Design Screens




## Visual Overview:

Clean, modern, and minimalist design, Soft gradient background (teal & white), Focused on clarity, comfort, and professionalism.

## Key Screen Components:

- **Top Bar:** Profile Icon (top-left), Settings Icon (top-right).
- **Temperature / Humidity Toggle:** Highlighted tab shows current selection
- **Health Status Indicator:** Displays "Good State" in bold font with gradient effect
- **Temperature Gauge:** Circular design showing 36°C with dynamic ring
- **Unit Switch:** Double-arrow icon toggles °C / °F for flexibility
- **Bottom Navigation Bar:** Icons: Home (active), Connectivity, Main Control, Control Mode, Maps.

# Light Mode




# Design Screens

A back button (arrow icon) is positioned on the top left, allowing users to return to the previous screen.

The left section contains a joystick-like circular control, likely used for analog movement control.

The right section features a directional button layout, which includes six directional arrows, indicating precise movement controls for the wheelchair.

# Light Mode Design Screens




The user profile screen shows a placeholder profile picture, "User Name" (ID: 23547858), and four card-like buttons: Personal Information, Medical History, Reminder, and Results and Examinations.

The Personal Information screen contains fields for Full name, Date, Age, Gender, Marital status, Residential address, Phone number, and Occupation, each with a placeholder value.

The Medical History screen shows a list of entries with columns for Date, Illness, and a delete icon. One entry is visible: Date: 2025-02-03, Illness: edasdasdasdas.

# Light Mode Design Screens



## Person Screen

### Screen Component:

Patient Avatar & Name.

### Feature Buttons:

Patient Information.  
Reminder.

Patient ID.

Medical History.  
Results and Examinations.

# Light

# Mode Design Screens

This screenshot shows the 'Personal Information' screen in light mode. The header is 'Personal Information'. Below it are several input fields with placeholder text: 'Full name: .....', 'Date: .....', 'Age: .....', 'Gender: .....', 'Marital status: .....', 'Residential address: .....', 'Phone number: .....', and 'Occupation: .....'. At the bottom is a navigation bar with icons for Home, Back, Share, and More.

2:32

Personal Information

Full name: .....

Date: .....

Age: .....

Gender: .....

Marital status: .....

Residential address: .....

Phone number: .....

Occupation: .....

## Personal Information Screen

### Screen Component:

Full name.

Marital status.

Occupation.

Date.


Residential address.

Age.

Gender.

Phone number.

# Light Mode Design Screens



## Medical History Screen


Each medical record is inside a rounded card, maintaining a consistent visual style. The details include:

- Date of record.
- Illness name (bolded for emphasis).
- Description (in a lighter font).

Each card has two action icons:

- Edit (pencil icon) – Allows modifying the record.
- Delete (trash icon) – Enables removal of the record.

# Light Mode Design Screens



## Reminder Screen

**Each Reminder record is inside a rounded card, maintaining a consistent visual style. The details include:**


- Date and Time: Displays when the reminder is set.
- Note: Provides a brief description of the reminder.

**Each card has two action icons:**

- Edit (pencil icon) – Allows modifying the record.
- Delete (trash icon) – Enables removal of the record.

# Light

# Mode Design Screens



## Examinations Screen


**Each Examinations record is inside a rounded card, maintaining a consistent visual style. The details include:**

- **Date Field:** Used to enter the date of the medical examination.
- **Illness Field:** Specifies the illness being examined.
- **Description Field:** Allows users to add additional notes about the examination.
- **Upload Files to Google Drive Button:** Enables users to upload medical reports or test results directly to Google Drive, enhancing cloud storage and accessibility of medical records.

**Each card has two action icons:**

- **Edit (pencil icon)** – Allows modifying the record.
- **Delete (trash icon)** – Enables removal of the record.

# Light Mode Design Screens



## Setting Page:

### Page Components:



Dark Mode Option.

Back Button  
(Top Left Corner).




Logout Option.

Team Member Page.




# Dark Mode Design Screens






# Dark Mode Design Screens






# Dark Mode Design Screens



# Dark Mode Design Screens



# Dark Mode Design Screens



UI

UI Code

Database

# Code

Code  
Flutter


# Flutter

# Why Flutter ?



# Flutter

# Block Diagram



# Block Diagram



UI

UI Code

Database

# Database


# Database

## Firebase

SVG



# Plugins Used



## SVG



## FIREBASE AUTH

Handles user authentication with email/password, Google, and other sign-in providers.



## URL LAUNCHER

Enables the app to open external URLs (like websites, email, or phone dialers) from within the app.



## FIREBASE CORE

The core library that connects your Flutter app to Firebase services. Required for any Firebase functionality.



## SVG



## FIREBASE AUTH

Handles user authentication with email/password, Google, and other sign-in providers.



## CLOUD FIRESTORE

Connects to Firebase Cloud Firestore, a real-time cloud-based NoSQL database.



## FIREBASE CORE

The core library that connects your Flutter app to Firebase services. Required for any Firebase functionality.



**SVG**



## FIREBASE AUTH

Handles user authentication with email/password, Google, and other sign-in providers.



## CLOUD FIRESTORE

Connects to Firebase Cloud Firestore, a real-time cloud-based NoSQL database.



## GOOGLE SIGN IN

Simplifies integrating Google Sign-In for authenticating users.



## FIREBASE CORE


The core library that connects your Flutter app to Firebase services. Required for any Firebase functionality.



# MOBILE & EMBEDDED LINK



# Block Diagram





# ESP 8266 Block Diagram

The ESP8266 is a low-cost Wi-Fi microcontroller module developed by Espressif Systems, widely used in IoT applications. It can connect to Wi-Fi networks and run standalone programs, making it ideal for smart devices and remote sensing.


## Key Features:


- **Multiple GPIO Pins:** Easy to interface with sensors, motors, and other hardware.
- **Arduino IDE Support:** Simple to program using Arduino libraries and tools.
- **Networking Protocols:** Supports TCP/IP, UDP, HTTP, MQTT for internet communication.





# Block Diagram






# Speech Recognition Model





# Speech Recognition Model

## Using Deep Learning

- We record the user's voice using a microphone .
- The audio is preprocessed and passed into a deep neural network to classify the spoken word.
- Once the command is recognized with prediction is longer than threshold ,it is sent to the STM32 microcontroller.



# STEPS FOR BUILDING THE MODEL



1

Data  
Collection

2

Data  
Preprocessing

3

Building  
Architecture  
model

4

Compilation  
process


5

Data  
training

6

Evaluation  
model

# STEPS FOR BUILDING THE MODEL



1

Data  
Collection


**Collect a high-quality dataset of audio samples by:**

- 1- sourcing publicly available speech datasets
- 2- creating your own recordings


**The labels are "forward , backward , left , right , and stop ". The audio files are in WAV format.**

<https://www.kaggle.com/competitions/tensorflow-speech-recognition-challenge/data>

# STEPS FOR BUILDING THE MODEL




**1**  
**Data  
Collection**



**2**  
**Data  
Preprocessing**




**3**  
**Building  
Architecture  
model**



**4**  
**Compilation  
process**




**5**  
**Data  
training**



**6**  
**Evaluation  
model**

# STEPS FOR BUILDING THE MODEL




2

Data  
Preprocessing


- 1) Sampling with Librosa Pad or truncate audio Exactly 1 second (8000 samples).
- 2) Data Augmentation: ( Add Noise - Time Shift - Time Stretch - Change Volume - Pitch Shift ).
- 3) Mel-Spectrogram Extraction:  
( 64-band mel-spectrogram- log-scaled dB spectrogram )
- 4) Normalization:  
( zero mean - unit variance - accelerates training ).
- 5) Data Preparation:  
Features stored  
Labels stored → one-hot encoded (5 classes)

# STEPS FOR BUILDING THE MODEL



1

Data  
Collection




2

Data  
Preprocessing



3

Building  
Architecture  
model




4

Compilation  
process



5


Data  
training



6

Evaluation  
model

# STEPS FOR BUILDING THE MODEL







3

Building  
Architecture  
model

- 1) Input Layer.
- 2) 1d Convolutions.
- 3) Batch Normalization && Pooling && Dropout Layers.
- 4) LSTM Layer.
- 5) Dense Layers.

# STEPS FOR BUILDING THE MODEL



1

Data  
Collection

2

Data  
Preprocessing

3

Building  
Architecture  
model

4

Compilation  
process


5

Data  
training

6

Evaluation  
model

# STEPS FOR BUILDING THE MODEL



4

Compilation  
process

## 1) Loss Function:

Choose an appropriate loss function, such as categorical cross-entropy, to measure the model's performance during training.


## 2) Optimizer:

Select an optimizer, like Adam or RMSprop, to update the model's weights and minimize the loss function.

## 3) Metrics:


Define relevant metrics, such as accuracy, to track the model's performance on the validation and test sets.

# STEPS FOR BUILDING THE MODEL



1

Data  
Collection




2

Data  
Preprocessing



3

Building  
Architecture  
model




4

Compilation  
process



5

Data  
training




6

Evaluation  
model





# STEPS FOR BUILDING THE MODEL



5  
Data  
training



# STEPS FOR BUILDING THE MODEL



1

Data  
Collection

2

Data  
Preprocessing

3

Building  
Architecture  
model

4

Compilation  
process


5

Data  
training

6

Evaluation  
model

# STEPS FOR BUILDING THE MODEL





6


Evaluation  
model

```
Epoch 30/50
517/517 - 6s 11ms/step - accuracy: 0.9698 - loss: 0.0920 - val_accuracy: 0.9685 - val_loss: 0.0963 - learning_rate: 1.2500e-04
Epoch 31/50
517/517 - 12s 14ms/step - accuracy: 0.9657 - loss: 0.0962 - val_accuracy: 0.9695 - val_loss: 0.0995 - learning_rate: 1.2500e-04
Epoch 32/50
517/517 - 9s 13ms/step - accuracy: 0.9713 - loss: 0.0865 - val_accuracy: 0.9700 - val_loss: 0.0985 - learning_rate: 1.2500e-04
Epoch 33/50
517/517 - 6s 12ms/step - accuracy: 0.9714 - loss: 0.0878 - val_accuracy: 0.9700 - val_loss: 0.0959 - learning_rate: 6.2500e-05
Epoch 34/50
517/517 - 11s 14ms/step - accuracy: 0.9730 - loss: 0.0823 - val_accuracy: 0.9702 - val_loss: 0.0977 - learning_rate: 6.2500e-05
130/130 - 0s 3ms/step - accuracy: 0.9713 - loss: 0.0999
```


Final Accuracy: 97.07% | Loss: 0.0997




# TESTING



# Block Diagram





# FACE AND EYE DIRECTION TRACKING MODEL



# Face Direction Tracking

## Methodology:

- Determines whether the user is looking forward, left, or right.
- Uses MediaPipe Face Mesh to detect landmarks (nose, eyes).
- Compares the nose's horizontal position relative to the eyes to assign a direction.

# Face Direction Tracking

## LOGIC


- Nose centered → Forward
- Nose shifted left → Left
- Nose shifted right → Right



Left



Forward





Right

# Face Direction Tracking

## Technologies Used:


- MediaPipe, OpenCV, Raspberry Pi Camera, Serial Communication



# Eye Tracking for Gaze Direction

## Methodology:


- Detects iris position inside the eye using image processing.
- **Steps include:**
  - Capture real-time video
  - Extract eye region
  - Apply grayscale, blur, and threshold
  - Use contour detection to find iris center



# Face Direction Tracking

## LOGIC

- Iris left → Left
- Iris right → Right
- Iris up → Forward
- Iris down → Down
- Iris centered → Stop







**MOBICARE**  
Smart Wheelchair App

Thanks